



125 SHERWOOD AVENUE  
FARMINGDALE, NY 11735

Phone: 631-501-1452 Fax: 631-501-1453

1-888-COLD-PAK | [www.nortechlabs.com](http://www.nortechlabs.com)

# EXAM DISPOSABLE & BLUE EASY SLEEVE SPEC SHEET

page 1 of 2

## I. PRODUCT IDENTIFICATION

### POLY-BOND POLYPROPYLENE NONWOVENS

All weights and colors of untreated polypropylene fabric

### CHEMICAL AND COMMON NAME

Polypropylene - family olefins

Poly-Bond™

(CH<sub>2</sub>-CH-CH<sub>3</sub>)<sub>n</sub>

### APPEARANCE AND ODOR

Point bonded nonwoven fabric, white, blue, and various colors

No odor

## II. HAZARDOUS INGREDIENTS AND EXPOSURE LIMITS

As we interpret the OSHA Hazard Communication Standard, polypropylene nonwoven fabric is not a hazardous material. According to the Code of Federal Regulations Title 29, Chapter XVII, Part 1910, Subpart 7. Paragraph (b) (5), the OSHA Hazard Communications Program does not apply to "(IV) Articles".

Nonwoven fabric exhibits performance parameters which meet all requirements as specified in the definition of an article. Therefore, the Hazard Communication Program does not apply to nonwoven fabric as an article.<sup>1</sup>

## III. TYPICAL PHYSICAL & CHEMICAL CHARACTERISTICS

Boiling Point	N/A
Vapor Pressure (mmHg)	N/A
Vapor Density	N/A
Volatile (VOL) %	Less than 0.4
Solubility in Water	Negligible
Specific Gravity	.83 - .92
Evaporation Rate	N/A

<sup>1</sup> A more detailed disclosure will be provided to appropriate personnel as privileged information upon request, in case of need for a specific application



### Nortech Labs, Inc.

*Providing green and innovative quality products to the medical, healthcare and cold chain shipping industry for over 50 years.*





125 SHERWOOD AVENUE  
FARMINGDALE, NY 11735

Phone: 631-501-1452 Fax: 631-501-1453

1-888-COLD-PAK | [www.nortechlabs.com](http://www.nortechlabs.com)

# EXAM DISPOSABLE & BLUE EASY SLEEVE SPEC SHEET

page 2 of 2

## Introducing The New “Deluxe SMS Apparel Material”

SMS Apparel is constructed from a unique Tri-Laminate Material Construction incorporating Spunbonded outer layers sandwiching an inner layer of Meltblown Fibers. This design construction delivers a fabric with High Tensile Strength and toughness that is also soft and drapeable, and in addition, provides a High Fluid Barrier that is Breathable. Now you can purchase a better material, design, construction and comfort at an affordable price.

### Fabric Features And Benefits:

- Breathable Barrier Protection With Comfort
- Heavy Duty Fabric Durability
- Multi-Layer Construction Blocks Liquid And Tiny Airborne Particles



### Nortech Labs, Inc.

*Providing green and innovative quality products to the medical, healthcare and cold chain shipping industry for over 50 years.*

